



A Novel Total-Body PET Scanner Using Xenon-Doped Liquid Argon Scintillator for Outstanding Detection

An application in medical physics of the DarkSide collaboration

Azam Zabihi, ASTROCENT/CAMK-PAN, Warsaw, Poland
On behalf of the 3Dπ TB-TOF-PET Collaboration

ASTROCENT



This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 952480



Republic of Poland




European Union
European Regional
Development Fund




LIDINE 2023 - Madrid

PET Principle



PET Annihilation Illustration



Three Types of Coincident Events

Liquid Xenon vs. Liquid Argon

And Benefit of Cryogenic

0.5%


Scintillator:	LAr	LXe	LAr + Xe	LYSO
Decay F/S (ns):	7/1600	4.3/22	~6/100*	41
Wavelength (nm):	128	175	~175**	420
Density (g/cm ³):	1.40	2.94	~1.40	7.1
Temperature (K):	87	162	87***	298
Photons/keV:	40	42	~41	28.5
Cost (US\$/kg):	~2	~2000	~2	~4

*Shorter slow decay time than the pure liquid argon.

**Scintillation light at a wavelength of 175 nm; Xe operates as a wavelength shifter (WLS).

***Operating at temperatures near the boiling point of argon eliminates the need for cooling and results in lower Dark Count Rate (DCR).

SiPM Dark Count Rate (DCR) vs. Temperature



<https://oar.princeton.edu/rt4ds/file/1663/1610.01915v1.pdf>

Reduction in the dark count rate (DCR), improves the timing capability of the devices and Signal-to-Noise Ratio (SNR)

3D π Overview

A **Total-body (TB)**, **Time of Flight (TOF)** PET scanner

- Xenon-doped Liquid Argon instead of Crystal scintillators
- Multiple detection layers
- Using Silicon Photomultipliers (SiPM)
- Double sided SiPM on scintillation

Geometry:


- 9 annulus detection layers
- Each layer has the scintillator sandwiched between two layers of SiPMs
- Each detection layer has ~18 mm LAr thickness
- PTFE supporting structure
- 2 m in length
- Geant4 simulations

Two configurations:

- ❖ **LAr+Xe, as the main focus**
- ❖ LAr+TPB (TetraPhenylButadiene: an organic WLS), for comparison

Geant4 Geometry Parameters

Parameter	Value
Inner radius (cm)	45
Outer radius (cm)	64
Length/AFOV (cm)	200
LAr thickness (cm)	16.2
Number of LAr layers	9
SiPM size (mm x mm)	10 x 10
Number of SiPMs	$\sim 1 \times 10^6$
Cryostat Thickness (mm)	6



3D π Geometry rendered in Fusion 360

Fluctuation Sources in Parallel and Perpendicular to LOR

Sources of fluctuation for construction of events **parallel to the LOR:**


1- Timing fluctuations:

- The statistical nature of scintillation photon detection,
- The time resolution of the SiPMs and their associated electronics.

2-Uncertainty resulting from imperfect scatter depth correction within the annular cylinder.

Sources of fluctuation for construction of events **perpendicular to the LOR:**

- The lattice size of the SiPM layout imposes a limitation on the resolution in each of the two directions perpendicular to the LOR.



Desirable timing resolution: $\sigma_{\text{SiPM}} < 60 \text{ ps}$.
Improved PDE (>40%) ensures acceptable resolution, independent of PDE.

In this work, we assume a σ_{SiPM} **60 ps** and the **PDE** is based on **ref [1]** at an over-voltage of 4 V.

National Electrical Manufacturers Association




NU 2-2018

A guide to characterize PET performance

Spatial Resolution of Single Point Sources

The spatial resolution of a system represents its ability to distinguish between two points after image reconstruction.




$$\theta_c = \text{atan}(y_c/x_c)$$

$$\rho = x \cos(\theta_c) + y \sin(\theta_c)$$


$$\alpha = z$$

$$\tau = y \cos(\theta_c) - x \sin(\theta_c)$$

Definition of radial, tangential, axial





Radial offsets



Spatial Resolution of Single Point Sources

Raw data (before image reconstruction stage)




LAr+Xe improves spatial resolution by more than a factor of 2.

Spatial Resolution of Single Point Sources

after image reconstruction: Filtered backprojection

Central Phantoms:	<u>1 cm radial position</u>			<u>10 cm radial position</u>			<u>20 cm radial position</u>		
<u>Scanner</u>	Radial [mm]	Tangent [mm]	Axial [mm]	Radial [mm]	Tangent [mm]	Axial [mm]	Radial [mm]	Tangent [mm]	Axial [mm]
3D π (LAR+Xe)	4.8	4.9	4.5	4.9	4.9	4.5	4.8	5.0	4.6
uExplorer (LYSO)	3.0	3.0	2.8	3.4	3.1	3.2	4.7	4.0	3.2



3D π is able to produce consistent and accurate images regardless of the location of the source.

Sensitivity


The sensitivity test measures the counts per second that the scanner measures for every unit of activity present in a source.

$$S_{\text{tot}} = \frac{R_{\text{CORR},0}}{A_{\text{cal}}}$$

S_{tot} : System Sensitivity

$R_{\text{CORR},0}$: The true coincidences count rate with no attenuation

A_{cal} : Line source radioactivity



Line source radial position in transaxial field of view	3Dπ LAr+Xe (200 cm AFOV) [kcps/MBq]	GE SIGNA PET/MR (25 cm AFOV) [kcps/MBq]	uExplorer PET/CT (192 cm AFOV) [kcps/MBq]
Center	564.0	21.8	174.0
10 cm radial offset	501.1	21.2	177.0

A higher system sensitivity indicates that the scanner can detect a larger fraction of the emitted photons, which allows for shorter scan times or lower radiotracer doses.

Noise Equivalent Count Rate (NECR)


$$NECR = \frac{T^2}{T + R + S}$$

Noise Equivalent Count Rate: ability to detect and accurately quantify true coincident counts while minimizing the impact of noise, (random, and scatter events.)




Source Distribution:

A solid right circular high density polyethylene cylinder with a line source.



Three types of coincident events

Noise Equivalent Count Rate (NECR)




---Background: Activity concentration used as the background in the Image Quality test
 ---Signal: Activity concentration used as the signal in the Image Quality test

$$NECR = \frac{T^2}{T + R + S}$$

Higher NECR with low activity indicates the possibility to reduce radioactive dose significantly

Image Quality

Measure **image contrast** and **background variability** using “hot” spheres in a uniform background.




Body (Cyan) and Test (Gold) phantom geometries
Activity concentration: 5.3 kBq/mL




Spheres (Red) phantom geometries
Activity concentration: $4 \times 5.3 \text{ kBq/mL} = 21.2 \text{ kBq/mL}$

Image Quality

3Dπ
scanned for 35 s.



uEXPLORER
scanned for 30 min.








Percent contrast=

$$\frac{\left(\frac{C_{H,j}}{C_{B,j}} \right) - 1}{\left(\frac{a_H}{a_B} \right) - 1} \times 100\%$$

$C_{H,j}$: the counts in the ROI for sphere j
 $C_{B,j}$: The average of the background ROI counts for sphere j
 a_B : The activity concentration in the background
 a_H : The activity concentration in the hot sphere

Comparison of NEMA Test Results

Scanner	Peak NECR [Mcps]	Activity concentration at peak [kBq/mL]	Sensitivity [kcps/MBq]	TOF resolution [ps]
 3Dπ (MC) (Preliminary)	~3	17.3*	560	163
	~3.5	30**		
 uEXPLORER TB-PET/CT	~1.5	17.3	174	412
 J-PET-TB (MC)	0.63	30	38	500
 GE SIGNA PET/CT	0.22	20.8	21.8	386
 VRain PET	0.14	9.8	25	229

The preliminary results demonstrate that our scanner system performance is comparable to commercial scanners.

*Activity concentration at peak NECR, uEXPLORER

**Activity concentration at peak NECR, J-PET

Ongoing activities: Cryogenics and Hardware

- Developing a hardware prototype of a LAr+Xe PET ring at Cagliari, Sardinia
 - ❖ A dedicated laboratory is being set up in Sardinia for conducting tests and commissioning the cryostat.
- INFN Torino developing a front end ASIC board


Next steps

Simulation/Software

- Experiment with newer reconstruction algorithms
- Improve Geant4 simulation
- Optimize detection layer geometry

Hardware

- Testing ALCOR board with LAr+Xe
- Ensuring stability and homogeneity of LAr+Xe
- Set up cryogenic infrastructure at Sardinia
- Develop the PET scanner prototype





Thank you

azabihi@camk.edu.pl